Instruction for Use

Ambu[®] aScope[™] 4 Broncho Family

Ambu[®] aScope[™] 4 Broncho Slim











Ambu[®] aScope[™] 4 Broncho Regular

Ambu[®] aScope[™] 4 Broncho Large

For use by trained clinicians/physicians only.
For in-hospital use.
For use with Ambu[®] aView[™] monitor.

CE 0086
Pat. Pending

Ambu is a registered trademark and aScope and aView are trademarks of Ambu A/S.

Contents	Page
English (Instruction for Use).....	4-13
Български (Указания за ползване).....	14-25
Česky (Návod k použití).....	26-35
Dansk (Brugsanvisning).....	36-45
Deutsch (Bedienungsanleitung).....	46-57
Ελληνικά (δηγίεις Χρήσεως).....	58-69
Español (Manual de instrucciones).....	70-79
Eesti (Kasutusjuhised).....	80-89
Suomi (Käyttöohje).....	90-99
Français (Mode d'emploi).....	100-109
Hrvatski (Upute za uporabu).....	110-119
Magyar (Használati útmutató).....	120-129
Italiano (Manuale d'uso).....	130-139
日本語 (使用法).....	140-149
Lietuviškai (Naudojimo instrukcijos).....	150-159
Latviski (Lietošanas instrukcija).....	160-169
Nederlands (Gebruiksaanwijzing).....	170-179
Norsk (Brukerveiledning).....	180-189
Polski (Instrukcja obsługi).....	190-201
Português (Manual de instruções).....	202-211
Português–Brazil (Manual de instruções).....	212-221
Română (Instrucțiuni de utilizare).....	222-231
Русский (Directions for use).....	232-243
Slovenčina (Návod na použitie).....	244-253
Slovenščina (Navodila za uporabo).....	254-263
Svenska (Instruktionshandbok).....	264-273
Türkçe (Kullanım talimatları).....	274-283
中文 (使用说明).....	284-291

1. Important Information – Read Before Use

Read these safety instructions carefully before using the aScope 4 Broncho. The Instruction for Use may be updated without further notice. Copies of the current version are available upon request. Please be aware that these instructions do not explain or discuss clinical procedures. They describe only the basic operation and precautions related to the operation of the aScope 4 Broncho. Before initial use of the aScope 4 Broncho, it is essential for operators to have received sufficient training in clinical endoscopic techniques and to be familiar with the intended use, warnings, cautions, indications and contraindications mentioned in these instructions. There is no warranty on the aScope 4 Broncho.

In this document *aScope 4 Broncho* refers to instructions which applies to the scope only and *aScope 4 Broncho system* refers to information relevant for the aScope 4 Broncho and the aView monitor and accessories. Unless specified otherwise, aScope 4 Broncho refers to all aScope 4 Broncho variants.

1.1. Intended Use

The aScope 4 Broncho endoscopes have been designed to be used with the aView monitor, endotherapy accessories and other ancillary equipment for endoscopy within the airways and tracheobronchial tree.

1.2. Indications for Use

The aScope 4 Broncho is for use in a hospital environment. The aScope 4 Broncho is a single-use device designed for use in adults. It has been evaluated for the following endotracheal tubes (ETT), double lumen tubes (DLT) and endoscopic accessories (EA) sizes:

	Minimum ETT inner diameter	Minimum DLT size	EA minimum working channel width
aScope 4 Broncho Slim 3.8/1.2	5.0 mm	35 Fr	up to 1.2 mm
aScope 4 Broncho Regular 5.0/2.2	6.0 mm	41 Fr	up to 2.0 mm
aScope 4 Broncho Large 5.8/2.8	7.0 mm	-	up to 2.6 mm

There is no guarantee that instruments selected solely using this minimum instrument channel width will be compatible in combination.

1.3. Warnings and Cautions

Failure to observe these warnings and cautions may result in patient injury or damage to the equipment. **Ambu is not responsible for any damage to the system or patient injury resulting from incorrect use.**

WARNINGS

1. aScope 4 Broncho is a single-use device and must be handled in a manner consistent with accepted medical practice for such devices in order to avoid contamination of the aScope 4 Broncho prior to insertion.
2. aScope 4 Broncho images must not be used as an independent diagnostic of any pathology. Physicians must interpret and substantiate any finding by other means and in the light of the patient's clinical characteristics.
3. Do not use active endoscopic accessories such as laser probes and electrosurgical equipment in conjunction with the aScope 4 Broncho, as this may result in patient injury or damage to the aScope 4 Broncho.
4. Do not use the aScope 4 Broncho if it is damaged in any way or if any part of the functional check (see section 4.1) fails.
5. Do not attempt to clean and reuse the aScope 4 Broncho as it is a single-use device. Reuse of the product can cause contamination, leading to infections.
6. The aScope 4 Broncho is not to be used when delivering highly flammable anaesthetic gases to the patient. This could potentially cause patient injury.
7. The aScope 4 Broncho is not to be used in a MRI environment.
8. Do not use the aScope 4 Broncho during defibrillation.
9. Only to be used by skilled physicians trained in clinical endoscopic techniques and procedures.

10. Do not use excessive force when advancing, operating or withdrawing the aScope 4 Broncho.
11. Patients should be adequately monitored at all times during use.
12. Always watch the live endoscopic image on the aView when advancing or withdrawing the aScope 4 Broncho, operating the bending section or suctioning. Failure to do so may harm the patient.
13. Do not use aScope 4 Broncho if the product sterilisation barrier or its packaging is damaged.
14. The distal end of the endoscope may get warm due to heating from the light emission part. Avoid long periods of contact between the tip of the device and the mucosal membrane as long, sustained contact with the mucosal membrane may cause mucosal injury.
15. Always make sure that any tube connected to the suction connector is connected to a suction device.
16. When withdrawing the endoscope, the distal tip must be in neutral and non-deflected position. Do not operate the bending lever, as this may result in injury to the patient and/or damage to the aScope 4 Broncho.
17. Do not advance or withdraw aScope 4 Broncho, or operate the bending section, while endoscopic accessories are protruding from the distal end of the working channel, as this may result in injury to the patient.
18. Always make sure that the bending section is in a straight position when inserting or withdrawing an endoscopic accessory in the working channel. Do not operate the bending lever and never use excessive force, as this may result in injury to the patient and/or damage to the aScope 4 Broncho.
19. Always perform a visual check according to the instructions in this Instructions for Use before placing the aScope 4 Broncho in a waste container.
20. Electronic equipment and the aScope 4 Broncho system may affect the normal function of each other. If the aScope 4 Broncho system is used adjacent to or stacked with other equipment, observe and verify normal operation of both the aScope 4 Broncho system and the other electronic equipment prior to using it. It may be necessary to adopt procedures for mitigation, such as reorientation or relocation of the equipment or shielding of the room in which it is used.
21. The aScope 4 Broncho consists of parts supplied by Ambu. These parts must only be replaced by Ambu authorised parts. Failure to comply with this may result in patient injury.
22. Be careful to check whether the image on the screen is a live image or a recorded image and verify that the orientation of the image is as expected.
23. To avoid risk of electric shock, the aScope 4 Broncho system must only be connected to a supply mains with protective earth. To disconnect the aScope 4 Broncho system from mains remove the mains plug from the wall outlet.
24. Always check compatibility with endotracheal tubes and double lumen tubes.
25. If any malfunction should occur during the endoscopic procedure, stop the procedure immediately and withdraw the endoscope.

CAUTIONS

1. Have a suitable backup system readily available for immediate use so the procedure can be continued if a malfunction should occur.
2. Be careful not to damage the insertion cord or distal tip when using sharp devices such as needles in combination with the aScope 4 Broncho.
3. Be careful when handling the distal tip of the insertion cord and do not allow it to strike other objects, as this may result in damage to the equipment. The lens surface of the distal tip is fragile and visual distortion may occur.
4. Do not exert excessive force on the bending section as this may result in damage to the equipment. Examples of inappropriate handling of the bending section include:
 - Manual twisting.
 - Operating it inside an ETT or in any other case where resistance is felt.
 - Inserting it into a preshaped tube or a tracheostomy tube with the bending direction not aligned with the curve of the tube.
5. US federal law restricts these devices for sale only by, or on the order of, a physician.
6. Keep the aScope 4 Broncho handle dry during preparation, use and storage.
7. Do not use a knife or other sharp instrument to open the pouch or cardboard box.
8. Secure the tubing properly on the suction connector before suction is applied.
9. If needed remove secretion or blood from the airway before and during the procedure. The suction function of any appropriate suction device can be used for this purpose.
10. Apply a vacuum of 85 kPa (638 mmHg) or less when suctioning. Applying too large a vacuum may make it difficult to terminate suctioning.


2. System Description

The aScope 4 Broncho can be connected to the aView monitor. For information about the aView monitor, please refer to the aView Instruction for Use.

2.1. System Parts

Ambu® aScope™ 4 Broncho
- Single use device:

Part numbers:

478001000 aScope 4 Broncho Large 5.8/2.8
477001000 aScope 4 Broncho Regular 5.0/2.2
476001000 aScope 4 Broncho Slim 3.8/1.2

aScope 4 Broncho Slim, aScope 4 Broncho Regular and aScope 4 Broncho Large are not available in all countries. Please contact your local sales office.

Product Name	Colour	Outer Diameter [mm]	Inner Diameter [mm]
aScope 4 Broncho Large 5.8/2.8	Orange	min 5.8; max 6.3	min 2.6
aScope 4 Broncho Regular 5.0/2.2	Green	min 5.0; max 5.5	min 2.0
aScope 4 Broncho Slim 3.8/1.2	Gray	min 3.8; max 4.2	min 1.2

Compatible devices: Ambu® aView™

Part numbers:

(Reusable)


405001000 Model no. JANUS-W08-R10
(SW versions v1.XX)*
405002000 Model no. JANUS2-W08-R10
(SW versions v2.XX)

For aView model no., please check the backside label on aView.

aView is not available in all countries. Please contact your local sales office.

* aScope 4 Broncho has limited image quality performance on aView model no. JANUS-W08-R10.


2.2. aScope 4 Broncho Parts

No.	Part	Function	Material
1	Handle	Suitable for left and right hand	MABS
2	Control lever	Moves the distal tip up or down in a single plane	POM
3	Working channel port	Allows for instillation of fluids and insertion of endoscopic accessories	MABS + Silicone
-	Working channel	Can be used for instillation of fluids, suction and insertion of endoscopic accessories	PU
4	Suction connector	Allows for connection of suction tubing	MABS
5	Suction button	Activates suction when pressed down	MABS
6	Tube connection	Allows for fixation of tubes with standard connector during procedure	TPE alloy
7	Insertion cord	Flexible airway insertion cord	PU
8	Bending section	Manoeuvrable part	PU
9	Distal end	Contains the camera, light source (two LEDs), as well as the working channel exit	Epoxy
7-8-9	Insertion portion	The ensemble of insertion cord, bending section and distal end	See above
10	Connector on aScope 4 Broncho cable	Connects to blue socket on aView	PVC
11	aScope 4 Broncho cable	Transmits the image signal to aView	PVC
12	Handle protection	Protects the suction connector during transport and storage. Remove before use.	PP
13	Protection pipe	Protects the insertion cord during transport and storage. Remove before use.	PP
14	Introducer	To facilitate introduction of Luer Lock syringes and soft endoscopic accessories through the working channel	PC
-	Packaging	Sterile barrier	Cardboard, Tyvek

Abbreviations: MABS (Methyl Acrylonitrile Butadiene Styrene), PU (Polyurethane), TPE (Thermoplastic Elastomer), PP (polypropylene), PC (Polycarbonate), POM (Polyoxymethylen).

3. Explanation of Symbols Used

Symbols for the aScope 4 Broncho devices	Indication
	Working length of the aScope 4 Broncho insertion cord.
 Max OD	Maximum insertion portion width (Maximum outer diameter).
 Min ID	Minimum working channel width (Minimum inner diameter).
 85°	Field of view.
	Do not use if the product sterilisation barrier or its packaging is damaged.
	This product is not made with natural rubber latex.
	Temperature limitation: between 10 °C (50 °F) and 40 °C (104 °F) in operating environment.
	Humidity limitation: relative humidity between 30 and 85% in operating environment.
	Atmospheric pressure limitation: between 80 and 109 kPa in operating environment.
	Manufacturer.
	Consult Instruction for Use.
	CE mark. The product complies with the EU Council directive concerning Medical Devices 93/42/EEC.
	Electrical Safety Type BF Applied Part.
	Use By, followed by YYYY-MM-DD.
	Sterile Product, Sterilisation by ETO.
	Single use product, do not reuse.
	Reference Number.
	Lot Number, Batch Code.
	UL Recognized Component Mark for Canada and the United States.
	Warning.

4. Use of aScope 4 Broncho

Numbers in gray circles below refer to illustrations on page 2.

4.1. Preparation and Inspection of aScope 4 Broncho

Visual inspection of the aScope 4 Broncho 1

1. Check that the pouch seal is intact. 1a
2. Make sure to remove the protective elements from the handle and from the insertion cord. 1b
3. Check that there are no impurities or damage on the product such as rough surfaces, sharp edges or protrusions which may harm the patient. 1c

Refer to the aView Instruction for Use for preparation and inspection of the aView monitor 2

Inspection of the Image

1. Connect aScope 4 Broncho to the aView by plugging the white aScope 4 Broncho connector with blue arrow into the corresponding blue female connector on the aView. Please be careful to align the arrows when inserting the aScope 4 Broncho connector into the aView monitor. 3
2. Verify that a live video image appears on the screen by pointing the distal end of aScope 4 Broncho towards an object, e.g. the palm of your hand. 4
3. Adjust the image preferences on the aView if necessary (please refer to the aView Instruction for Use).
4. If the object cannot be seen clearly, wipe the lens at the distal end using a sterile cloth.

Preparation of aScope 4 Broncho

1. Carefully slide the bending control lever forwards and backwards to bend the bending section as much as possible. Then slide the bending lever slowly to its neutral position. Confirm that the bending section functions smoothly and correctly and returns to a neutral position. 5a
2. Using a syringe insert 2ml of sterile water into the working channel port (if applying a Luer Lock syringe use the enclosed introducer). Press the plunger, ensure that there are no leaks, and that water is emitted from the distal end. 5b
3. If applicable, prepare the suction equipment according to the supplier's manual. 5c
Connect the suctioning tube to the suction connector and press the suction button to check that suction is applied.
4. If applicable, verify that endoscopic accessory of appropriate size can be passed through the working channel without resistance. The enclosed introducer can be used to facilitate the insertion of soft accessories such as microbiology brushes. 5d
5. If applicable, verify that endotracheal tubes and double lumen tubes are compatible with aScope 4 Broncho before starting the procedure.

4.2. Operating the aScope 4 Broncho

Holding the aScope 4 Broncho and manipulating the tip 6a

The handle of the aScope 4 Broncho can be held in either hand. The hand that is not holding the aScope 4 Broncho can be used to advance the insertion cord into the patient's mouth or nose. Use the thumb to move the control lever and the index finger to operate the suction button. The control lever is used to flex and extend the tip of the aScope 4 Broncho in the vertical plan. Moving the control lever downward will make the tip bend anteriorly (flexion). Moving it upward will make the tip bend posteriorly (extension). The insertion cord should be held as straight as possible at all times in order to secure an optimal tip bending angle.

Tube connection 6b

The tube connection can be used to mount ETT and DLT with an ISO connector during intubation.

Insertion of the aScope 4 Broncho 7a

Lubricate the insertion cord with a medical grade lubricant when the aScope 4 Broncho is inserted into the patient. If the camera image of the aScope 4 Broncho becomes unclear the tip can be cleaned by gently rubbing the tip against the mucosal wall or remove the aScope 4 Broncho and clean the tip. When inserting the aScope 4 Broncho orally, it is recommended to use a mouthpiece to protect the aScope 4 Broncho from being damaged.

Instillation of fluids 7b

Fluids can be instilled through the working channel by inserting a syringe into the working channel port at the top of the aScope 4 Broncho. When using a Luer Lock syringe, use the included introducer. Insert the syringe completely into the working channel port or the introducer and press the plunger to instill fluid. Make sure you do not apply suction during this process, as this will direct the instilled fluids into the suction collection system. To ensure that all fluid has left the channel, flush the channel with 2ml of air.

Aspiration 7c

When a suction system is connected to the suction connector, suction can be applied by pressing the suction button with the index finger. If the introducer and/or an endoscopic accessory is placed inside the working channel note that the suction capability will be reduced. For optimal suction capability it is recommended to remove the introducer or syringe entirely during suction.

Insertion of endoscopic accessories 7d

Always make sure to select the correct size endoscopic accessory for the aScope 4 Broncho (See section 1.2). Inspect the endoscopic accessory before using it. If there is any irregularity in its operation or external appearance, replace it. Insert the endoscopic accessory into the working channel port and advance it carefully through the working channel until it can be seen on the aView. The enclosed introducer can be used to facilitate the insertion of soft accessories such as microbiology brushes.

Withdrawal of the aScope 4 Broncho 8

When withdrawing the aScope 4 Broncho, make sure that the control lever is in the neutral position. Slowly withdraw the aScope 4 Broncho while watching the live image on aView.

4.3. After Use

Visual check 9

1. Are there any missing parts on the bending section, lens, or insertion cord? If yes, then take corrective action to locate the missing part(s).
2. Is there any evidence of damage on the bending section, lens, or insertion cord? If yes, then examine the integrity of the product and conclude if there are any missing parts.
3. Are there cuts, holes, sagging, swelling or other irregularities on the bending section, lens, or insertion cord? If yes, then examine the product to conclude if there are any missing parts.

In case of corrective actions needed (step 1 to 3) act according to local hospital procedures. The elements of the insertion cord are radio opaque.

Final steps

Disconnect the aScope 4 Broncho from the aView 10. The aScope 4 Broncho is a single use device. The aScope 4 Broncho is considered infected after use and must be disposed of in accordance with local guidelines for collection of infected medical devices with electronic components. Do not soak, rinse, or sterilize this device as these procedures may leave harmful residues or cause malfunction of the device. The design and material used are not compatible with conventional cleaning and sterilization procedures.

5. Technical Product Specifications

5.1. Standards Applied

The aScope 4 Broncho function conforms with:

- Council Directive 93/42/EEC concerning Medical Devices.
- IEC 60601-1 ed 2 Medical electrical equipment – Part 1: General requirements for safety.
- IEC 60601-1 ed 3.1 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance.
- IEC 60601-2-18 Medical electrical equipment – Part 2-18: Particular requirements for the safety of endoscopic equipment.
- ISO 8600-1: Optics and photonics - Medical endoscopes and endotherapy devices - Part 1: General requirements.
- IEC 60601-1-2: Medical electrical equipment – Part 1-2 General requirements for safety – Collateral standard: Electromagnetic compatibility - Requirements for test.
- ISO 10993-1: Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
- ISO 594-1: Conical fittings with 6% (Luer) taper for syringes, needles and certain other medical equipment – Part 1: General requirements.

5.2. aScope 4 Broncho Specifications

Insertion portion	aScope 4 Broncho Slim	aScope 4 Broncho Regular	aScope 4 Broncho Large
Bending section ¹ [°]	180 ↑, 180 ↓	180 ↑, 180 ↓	180 ↑, 160 ↓
Insertion cord diameter [mm, (")]	3.8 (0.15)	5.0 (0.20)	5.8 (0.23)
Distal end diameter [mm, (")]	4.2 (0.16)	5.4 (0.21)	6.2 (0.24)
Maximum diameter of insertion portion [mm, (")]	4.3 (0.17)	5.5 (0.22)	6.3 (0.25)
Minimum endotracheal tube size (ID) [mm]	5.0	6.0	7.0
Minimum double lumen tube size (ID) [Fr]	35	41	-
Working length [mm, (")]	600 (23.6)	600 (23.6)	600 (23.6)
Channel	aScope 4 Broncho Slim	aScope 4 Broncho Regular	aScope 4 Broncho Large
Minimum instrument channel width ² [mm, (")]	1.2 (0.047)	2.0 (0.079)	2.6 (0.102)
Storage and transportation	aScope 4 Broncho Slim/Regular/Large		
Transportation temperature [°C, (°F)]	10 ~ 40 (50 ~ 104)		
Recommended storage temperature ³ [°C, (°F)]	10 ~ 25 (50 ~ 77)		
Relative humidity [%]	30 ~ 85		
Atmospheric pressure [kPa]	80 ~ 109		

Optical System	aScope 4 Broncho Slim/Regular/Large
Field of View [°]	85
Depth of Field [mm]	6 - 50
Illumination method	LED
Suction connector	
Connecting tube ID [mm]	Ø7 +/- 1
Sterilisation	
Method of sterilisation	ETO
Operating environment	
Temperature [°C, (°F)]	10 ~ 40 (50 ~ 104)
Relative humidity [%]	30 ~ 85
Atmospheric pressure [kPa]	80 ~ 109
Altitude [m]	≤ 2000

1. Please be aware that the bending angle can be affected if the insertion cord is not kept straight.
2. There is no guarantee that accessories selected solely using this minimum instrument channel width will be compatible in combination.
3. Storage under higher temperatures may impact shelf life.

6. Trouble Shooting

If problems occur with the system, please use this trouble shooting guide to identify the cause and correct the error.

Problem	Possible cause	Recommended action
No live image on the left side of the screen but User Interface is present on the display or the image shown to the left is frozen	aScope 4 Broncho not connected to aView.	Connect an aScope 4 Broncho to the blue port on aView.
	aView and aScope 4 Broncho have communication problems.	Restart aView by pressing the power button for at least 2 seconds. When aView is off restart by pressing power button once more.
	aScope 4 Broncho is damaged.	Replace the aScope 4 Broncho with a new one.
	A recorded image is shown in the yellow file management tab.	Return to live image by pressing the blue live image tab or restart aView by pressing the power button for at least 2 seconds. When aView is off restart by pressing power button once more.
Low picture quality	Blood, saliva etc. on the lens (distal tip).	Gently rub the distal tip against the mucosa. If the lens cannot be cleaned this way remove the aScope 4 Broncho and wipe the lens with sterile gauze.
Absent or reduced suction capability or difficulty in inserting endoscopic accessory through the channel	Channel blocked.	Clean the working channel using a cleaning brush or flush the working channel with sterile saline using a syringe. Do not operate the suction valve when instilling fluids.
	Suction pump is not turned on or not connected.	Turn the pump on and check the suction line connection.
	Suction valve is damaged.	Prepare a new aScope 4 Broncho.
	Endoscopic accessory inserted in working channel (applicable if suction is absent or reduced).	Remove endoscopic accessory. Check that the accessory used is of the recommended size.
	Bending section not in neutral position.	Move bending section into neutral position.
	Soft endoscopic accessory difficult to pass through working channel seal.	Use the enclosed introducer.

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